SURGICAL ELEMENTS / SURGICAL CARE AUSTRALIA

ABN: 72 783 636 436 TGA. ARTG no: 225282

29 Exchange Ave, Harrisdale Western Australia 6112

info@surgicalcareaustralia.com.au info@surgicalelements.com

Ph: +61 (0) 8 9397 9988 M: +61 (0) 419 923 048 www.surgicalcareaustralia.com.au www.surgicalelements.com

Surgical Care Australia is a high quality "Surgical Patient Positioning Device Manufacturer" specialising in Orthopaedics. We have over 25 years design, development and manufacture experience in the Biomedical Engineering field. Surgical Care Australia works closely with teaching hospitals, surgeons, specialists and health professionals to fully appreciate and understand their needs and to ensure that our designs and products rigorously meet their requirements.

At Surgical Care Australia, we are committed to excellence in design, extremely high quality and to continued product innovation. We strive to deliver exceptional customer service and to constantly adapt to the changing needs of our profession, our clients and their patients.

To order, please contact Albert Hansma on 0419 923 048 or email info@surgicalcareaustralia.com.au *Please check our website as products are added regularly.*






TRC

Standard operation table

rail clamp

NOTE: specify your table make/model and rail size prior to ordering



AK-05-25 Dual strap 25mm Knee Brace Arthroscopy

•ACL

Knee Arthroscopy

Trauma

AK-05-50 Single strap 50mm

* trusted and durable design



AK-06-330 standard AK-06-420 bariatric Knee Brace Arthroscopy

·ACL

Knee Arthroscopy

Trauma

* only one person needed to set up knee/leg

* lightweight



HAT-300 300mm paediatric HAT-600 600mm standard


HAT-700 700mm tall patient, long arm Hand/arm surgical table

Radio translucent carbon fibre + vinyl mattress

Carbon fiber

Surgical hand / arm table

* the most rigid, stable board on the market



PC-90 90 degree 16mm round post clamp


Multi-purpose support post clamp

Multi-purpose




PK-01 Popliteal roll • Knee 80mm round roll • Trauma

* integral skin foam / impact absorption / anti-fluid absorption




LBK-02 Back support lateral • Hip (including bracket) • Shoulder • Multi-purpose


* integral skin foam / impact absorption / anti-fluid absorption



FR-01 Footrest with flat bar post • Total knee 10 x 20mm • ACL knee trauma

* integral skin foam / impact absorption / anti-fluid absorption










* Left: surgical patient positioning devices assembled on workstation






Function

LBK-01

Lateral back support


- Shoulder
- Patient safety barrier
- Trauma



SP-01

Straight stainless steel post

Multi-purpose



AP-01

Angled stainless steel post

Multi-purpose



* Left: straight and angled steel posts assembled on workstation






MCD-S (small) McDonald Splint (glove)

MCD-M (medium) MCD-L (large) Shoulder Arthroscopy

Function

Used in combination with:

- · SGA shoulder gantry
- Sand bags /counter weights 2.5/5kg



PSB-BC

Patient side board Atree Board + Cover Operating table safety

barrier

* unbreakable internal aluminium board



OWS-01

Orthopaedic Workstation Cart

* Overall Height: 950mm

Width: 550mm Length:700mm

* Basket draw H: 280mm

W: 420mm

L: 500mm

* Wheel diameter: 100mm

^{*} comfortable splint / machine washable



AE-01 Elbow arthroscopy brace

* no fuss - easy to use



SGA-01 Shoulder gantry Shoulder arthroscopy

Used in combination with:


- · MCD-S
- · MCD-M
- · MCD-L
- Counter weights / sand bag 2.5/5kg

* fast and easy set-up



Sandbags

* dual vinyl lining




CDP-01 Drape screen pole clamp Max. pole diameter 18mm



Foam Blocks - various

- · Lateral u-channel
- Head-box
- · Ankle/leg rest



PSP-01

- Plain stainless steel pole with 3 flat sides
- · Length 600mm

For ordering and further information, please contact Albert on 0419 923 048 or email info@surgicalcareaustralia.com.au

Check our website for product updates: www.surgicalcareaustralia.com.au

Warranty statement: All products are sold with a 12 month warranty, against faulty material and workmanship, effective from shipping date. A suspect product will be repaired or replaced free of charge, should a claim be deemed fair and supported by this supplier and client. Client must contact this supplier with the product model number, serial number and purchase date, for an immediate warranty investigation. Claims cannot be supported where it is clear that the product has been used in a way other than its intended purpose.